

## PARTIALLY EDENTULOUS CBCT SCAN OF BOTH ARCHES

### 1. Patient Positioning:

- i** - Breathe through your nose
- Bite gently on cotton rolls
- Close your eyes and place tongue at the back of the mouth

### 2. Vital anatomy to be captured:

- i** - Superior: Sinus (below orbit)
- Inferior: Chin and entire mandible



### 3. CBCT General Settings: Full 13cm-20cm, 0.3-0.4 Voxel

### 4. Export and Upload CBCT Scan

- i** - Export the CBCT scan in a multi-file DICOM format to the desktop
- Right click on the folder, go to "Send To", and select "Compressed (Zipped) Folder"
- Upload the zip / compressed folder to your Implant Concierge case

### 5. Ship the diagnostic model to Implant Concierge.

- i** - For more information on exporting DICOM and printing a shipping label, please visit:  
[www.implantconcierge.com/support](http://www.implantconcierge.com/support)